Flow Cell Wash Kit (EXP-WSH004)

Flow Cell Wash Kit (EXP-WSH004) Version: WFC_9120_v1_revI_08Dec2020 Last update: 12/04/2023 Flow Cell Number: DNA Samples:	Oxford NANOPORE Technologies
Flow Cell Number: DNA Samples:	NOTES/OBSERVATIONS
☐ Slowly add 500 µl of Storage Buffer (S) through the flow cell priming port.	
☐ Close the priming port.	
☐ Using a P1000, remove all fluid from the waste channel through Waste port 1. As both the flow cell priming port and SpotON sample port are closed, no fluid should leave the sensor array area.	
It is vital that the flow cell priming port and SpotON sample port are closed to prevent air from being drawn across the sensor array area, which would lead to a significant loss of sequencing channels.	
☐ The flow cell can now be stored at 4-8°C.	
When you wish to reuse the flow cell, remove the flow cell from storage, and allow it to warm to RT for \sim 5 minutes.	
After performing a flow cell wash, we recommend using the first pore scan to check number of available nanopores. After performing a flow cell wash, we recommend using the first pore scan to check number of available nanopores. Stop your sequencing run, prime your flow cell and load the library and start a new sequencing run OR pause your sequencing run, prime your flow cell and load the library and restart the sequencing run Wahing 참고 영상 https://www.google.com/search?q=flow+cell+washing&rlz =1C1GCEU_koKR1050KR1050&sxsrf=APwXEddgjBmjF_J5jca yQmQG6w7AfdB- EQ%3A1682921801521&ei=SVIPZOK4H8_d2roPv8iC- AU&oq=flow+cell+lwash&gs_lcp=Cgxnd3Mtd2l6LXNlcnAQA RgBMgclABANEIAEMgclABANEIAEMgclABANEIAEMgclABA NEIAEMgclABANEIAEMgclABANEIAEMgclABANEIAEMgclABA	
4QDTIGCAAQHhANMgYIABAeEA06CggAEEcQ1gQQsAM6BA gjECc6CAgAEIAEEMsBOgQIABAeOgYIABAeEA86BggAEAUQ HjoGCAAQCBAeOgUIIRCgAUoECEEYAFCEBFjdE2CMKmgBcA F4AIABiAGIAdsEkgEDMC41mAEAoAEByAEJwAEB&sclient=g ws-wiz- serp#fpstate=ive&vld=cid:db69ba0a,vid:KQBOypLcnOE	

nanoporetech.com Page 11/ 3