

Ligation sequencing amplicons - Native Barcoding Kit 24 V14 (SQK-NBD114.24)

Version: NBA_9168_v114_revE_15Sep2022
Last update: 06/01/2023

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><input type="checkbox"/> Keep the tube on the magnetic rack and wash the beads with 700 µl of freshly prepared 80% ethanol <u>without disturbing the pellet</u>. Remove the ethanol using a pipette and discard.</p> <p><input type="checkbox"/> Repeat the previous step.</p> <p><input type="checkbox"/> Spin down and place the tube back on the magnetic rack. Pipette off any residual ethanol. <u>Allow the pellet to dry for ~30 seconds, but do not dry the pellet to the point of cracking.</u></p> <p><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend the pellet in 35 µl Nuclease-free water by gently flicking.</p> <p><input type="checkbox"/> Incubate for 10 minutes at 37°C. Every 2 minutes, agitate the sample <u>by gently flicking</u> for 10 seconds to encourage DNA elution.</p> <p><input type="checkbox"/> Pellet the beads on a magnetic rack until the eluate is clear and colourless. Remove and retain 35 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.</p>	<p>37도 10분간 incubation, 2분마다 tapping</p> <p>다음 step에서 30ul 필요하므로 갈색 AXP가 딸려오지 않도록 주의</p>
<p>Quantify 1 µl of eluted sample using a Qubit fluorometer.</p>	<p>DNA loss가 생기는 지 매 step마다 농도 측정하여 확인, loss가 심하면 sample을 버리지 않고 현재 step을 다시해야함.</p>
<p>Take forward the barcoded DNA library to the adapter ligation and clean-up step. However, at this point it is also possible to store the sample at 4°C overnight.</p>	<p>실험 Stop point (2)</p>
<p>[3] Adapter ligation and clean-up</p>	<p>시퀀싱을 위한 어댑터를 붙여주는 과정 50min. 바코드 붙인 샘플을 동일한 농도로 pooling 후 진행. Adaptor ligation 후 최종 10-20 fmol이 나와야 하므로 loss가 안생기도록 주의하기. 농도가 낮은 경우 2step에서 다시 진행.</p>
<p>IMPORTANT</p> <p><input type="checkbox"/> The Native Adapter (NA) used in this kit and protocol <u>is not interchangeable with other sequencing adapters.</u></p>	
<p>Prepare the <u>NEBNext Quick Ligation Reaction Module</u> according to the manufacturer's instructions, and place on ice:</p> <p><input type="checkbox"/> Thaw the reagents at RT.</p> <p><input type="checkbox"/> Spin down the reagent tubes for 5 seconds.</p> <p><input type="checkbox"/> Ensure the reagents are fully mixed by performing 10 full volume pipette mixes.</p>	
<p>IMPORTANT</p> <p><input type="checkbox"/> Do not vortex the Quick T4 DNA Ligase.</p>	
<p><input type="checkbox"/> Spin down the <u>Native Adapter (NA)</u> and <u>Quick T4 DNA Ligase</u>, pipette mix and place on ice.</p> <p><input type="checkbox"/> Thaw the <u>Elution Buffer (EB)</u> at RT, mix by vortexing, spin down and place on ice.</p>	
<p>IMPORTANT</p> <p>Depending on the wash buffer (LFB or SFB) used, the clean-up step after adapter ligation is designed to either enrich for DNA fragments of >3 kb, or purify all fragments equally.</p> <p><input type="checkbox"/> To enrich for DNA fragments of 3 kb or longer, use Long Fragment Buffer (LFB)</p> <p><input type="checkbox"/> To retain DNA fragments of all sizes, use Short Fragment Buffer (SFB)</p>	<p>Metagenome의 경우 SFB 사용함.</p>
<p><input type="checkbox"/> Thaw either Long Fragment Buffer (LFB) or Short Fragment Buffer (SFB) at RT, mix by vortexing, spin down and place on ice.</p>	

2nd