

Flow Cell Wash Kit (EXP-WSH004)

Version: WFC_9120_v1_revI_08Dec2020
Last update: 12/04/2023



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 µl, and make sure that the array of pores are covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores. <input type="checkbox"/> Load 400 µl of the prepared Flow Cell Wash Mix into the flow cell priming port, avoiding the introduction of air. <input type="checkbox"/> Close the flow cell priming port and wait for 60 minutes. <input type="checkbox"/> Before removing the waste fluid a second time, ensure that the flow cell priming port cover and SpotON sample port cover are closed, as indicated in the figure below. <input type="checkbox"/> Using a P1000, remove all fluid from the waste channel through Waste port 1. As both the flow cell priming port and SpotON sample port are closed, no fluid should leave the sensor array area. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> It is vital that the flow cell priming port and SpotON sample port are closed to prevent air from being drawn across the sensor array area, which would lead to a significant loss of sequencing channels. 	
<p>Follow one of the two options described in the next steps of the protocol</p> <p>To run a second library on a MinION/GridION flow cell straight away</p>	<p>바로 사용할 경우 QC 체크 필요. 하지만 4도에 어느정도 보관 후 사용하는 것을 추천 (온도가 올라가서 pore에 무리가 올 수 있으므로)</p>
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> The buffers used in this process are incompatible with conducting a Flow Cell Check step prior to loading a subsequent library. However, the first pore scan after a sequencing run has started will report the number of nanopores available. <input type="checkbox"/> To run a second library straight away, follow the instructions in the "Priming and loading the flow cell" section of your library preparation protocol. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> When priming a flow cell after a nuclease wash with the Flow Cell Wash Kit, it is vital to wait five minutes between the priming mix flushes to ensure effective removal of the nuclease. 	
<p>To store the MinION/GridION flow cell for later use</p>	<p>4도에 보관할시 바로 이 step</p>
<p>Storage Buffer (S) can be used to flush flow cells for storage for later use or to check number of available nanopores before loading another library.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Thaw one tube of Storage Buffer (S) at RT. <input type="checkbox"/> Mix contents thoroughly by pipetting and spin down briefly. <input type="checkbox"/> Rotate the flow cell priming port cover clockwise so that the priming port is visible. <p>After opening the priming port, check for a small air bubble under the cover. Draw back a small volume to remove any bubbles:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Set a P1000 pipette to 200 µl. <input type="checkbox"/> Insert the tip into the flow cell priming port. <input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, or until you can see a small volume of buffer/liquid entering the pipette tip. <input type="checkbox"/> Visually check that there is continuous buffer from the flow cell priming port across the sensor array. 	<p>액체가 조금 나왔다 싶으면 stop !!! 기포가 들어가면 안됩니다.</p>