

Ligation sequencing amplicons - Native Barcoding Kit 24 V14 (SQK-NBD114.24)

Version: NBA_9168_v114_revE_15Sep2022
Last update: 06/01/2023



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
Remove and retain 10 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube. <input type="checkbox"/> Dispose of the pelleted beads	
1st Quantify 1 µl of each eluted sample using a Qubit fluorometer.	
Take forward an equimolar mass of samples to be barcoded and pooled forward into the native barcode ligation step. However, at this point it is also possible to store the sample at 4°C overnight .	실험 Stop point (1)
[2] Native barcode ligation	바코드를 붙여주는 과정 60min (NB01-24까지 있음. 한 flowcell에 중복되게 사용하지 않기)
<p>Prepare the NEB Blunt/TA Ligase Master Mix according to the manufacturer's instructions, and place on ice:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Thaw the reagents at RT. <input type="checkbox"/> Spin down the reagent tubes for 5 seconds. <input type="checkbox"/> Ensure the reagents are fully mixed by performing 10 full volume pipette mixes. <input type="checkbox"/> Thaw the EDTA at RT, mix by vortexing, spin down and place on ice. <input type="checkbox"/> Thaw the Native Barcodes (NB01-24) required for your number of samples at RT. Individually mix the barcodes by pipetting, spin down, and place them on ice. <input type="checkbox"/> Select a unique barcode for each sample to be run together on the same flow cell. Up to 24 samples can be barcoded and combined in one experiment. <p>In clean 0.2 ml PCR-tubes or a 96-well plate, add the reagents in the following order per well: total(20ul)</p> <ul style="list-style-type: none"> <input type="checkbox"/> 7.5 µl End-prepped DNA <input type="checkbox"/> 2.5 µl Native Barcode (NB01-24) <input type="checkbox"/> 10 µl Blunt/TA Ligase Master <p>Mix</p> <ul style="list-style-type: none"> <input type="checkbox"/> Ensure the reaction is thoroughly mixed by gently pipetting and spin down briefly. <input type="checkbox"/> Incubate for 20 minutes at RT. <input type="checkbox"/> Add 2 µl of EDTA to each well and mix thoroughly by pipetting and spin down briefly. <p><input type="checkbox"/> Pool the barcoded samples in a 1.5 ml Eppendorf DNA LoBind tube.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Resuspend the AMPure XP Beads (AXP) by vortexing. <input type="checkbox"/> Add AMPure XP Beads (AXP) to the pooled reaction, and mix by pipetting for a 0.4X clean. <input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 10 minutes at RT. <input type="checkbox"/> Prepare 2 ml of fresh 80% ethanol in Nuclease-free water. <p>Spin down the sample and pellet on a magnet for 5 minutes. Keep the plate on the magnetic rack until the eluate is clear and colourless, and pipette off the supernatant.</p>	<p>한 flowcell에 중복되게 사용하지</p> <p>바코드 하나에 F/R 동시에 들어가 있음. 한 샘플당 바코드 하나</p> <p>adding EDTA to stop the reaction</p> <p>바코딩할 샘플당 8ul AXP 첨가</p> <p>2nd</p> <p>홀라믹서가 없을 시 손으로 inverting</p>